


Zwitterion Chromatography – ZIC™

A novel technique, with unique selectivity,
suitable for preparative scale separations?

PhD Einar Pontén

What is Zwitterion Chromatography?



Our definition:

Liquid chromatography on a stationary phase having covalently bound zwitterionic functional groups.

The ZIC™ stationary phase

- ❖ Balanced stoichiometry and a **zero net charge**
- ❖ **Weak electrostatic interactions** with charged analytes
- ❖ Separation may be performed in totally **aqueous buffers**
- ❖ **High recovery** – absence of hydrophobic interactions
- ❖ Selectivity benefit by **charge** and **hydrophilicity**
- ❖ Suitable for *Hydrophilic Interaction Chromatography*
→ **ZIC™-HILIC Columns**

The ZIC™-HILIC stationary phase

- ❖ Permanent zwitterion in the practical pH range
- ❖ Stoichiometric molar balance sulfur:nitrogen (~1:1)
- ❖ Spherical porous silica particles (typically 5 µm, 200Å)
- ❖ Carbon content; ~10 %

Hydrophilic Interaction Chromatography (HILIC)


Isocratic or gradient elution using a mobile phase having an initial high content (>70%) of acetonitrile or low alcohols on a hydrophilic stationary phase, e.g., anion exchangers.

How do I use these columns?

The ZIC™ mode

- ❖ Separation by weak electrostatic interactions in aqueous eluents without organic solvents.
- ❖ Separation of proteins using gradient elution and typically lower electrolyte concentration than for traditional ion-exchange phases.
- ❖ Protein selectivity comparable to cation-exchange chromatography.

Amino acid retention in ZIC™ mode




- ❖ Low retention of AA with hydrophobic side groups (Leu)
- ❖ Low retention of hydrophilic AA without charged side groups (Gly, Met)
- ❖ Retention of basic AA due to electrostatic interaction with the zwitterionic group (Lys, His)

How do I use these columns?

The ZIC™-HILIC mode


- ❖ An initially high content of organic solvent (e.g., alcohol, acetonitrile) in the eluent. Gradient elution by increasing the content of water and eventually raising electrolyte concentration. Volatile buffers may be used.
- ❖ Unique selectivity from hydrophilic interactions superimposed onto weak electrostatic interactions.
- ❖ Suitable for peptides, carbohydrates, natural products, plant hormones, amino acids, and numerous more or less polar analytes.

Amino acid retention in ZIC™-HILIC mode



- ❖ Retention of AA increases with the concentration of acetonitrile in the eluent
- ❖ Retention increases with AA hydrophilicity
Gly > Met > Leu

ZIC™-HILIC separation of amino acids



Retention determined by:

- ❖ Hydrophilicity
- ❖ Charge

Linear gradient 4 %/min of Eluent B


Eluent A: 85 %(v/v) acetonitrile, 15 %(v/v)
10 mM KH₂PO₄ pH 4.5

Eluent B: 40 mM KH₂PO₄ pH 4.5

Flow rate: 0.8 mL/min, UV detection 214 nm
Injection volume: 20 µL (2.5 mg/mL each of
gly, leu and lys; 0.15 mg/mL His).

ZIC™-HILIC column (4.6×150 mm, 200Å, 5 µm)

ZIC™ separation of peptides



Isocratic separation by:

❖ Charge

Bradykinin:

Arg-Pro-Pro-Gly-Phe-Ser-Pro-Phe-Arg

Angiotensin II:

Asp-Arg-Val-Tyr-Ile-His-Pro-Phe


Eluent: formic acid at pH 2.6

Flow rate: 0.8 mL/min, UV detection 214 nm

Injection volume: 20 µL.

ZIC™-Si column (3 × 150 mm, 200Å, 10 µm).

Separation of peptides by ZIC™-HILIC



Retention determined by:

- ❖ Hydrophilicity
- ❖ Charge

Linear gradient 10-40 % of Eluent B in 25 minutes


Eluent A: 85 %(v/v) acetonitrile, 15 %(v/v) 10 mM KH₂PO₄ pH 4.5

Eluent B: 20 mM KH₂PO₄ pH 4.5

Flow rate: 0.8 mL/min, UV detection 214 nm


Injection volume: 20 µL

ZIC™-HILIC column (4.6×150 mm, 5 µm)



.....or by RPLC

Less successful due to peptide hydrophilicity and/or charge




**Linear gradient 0-60 % of Eluent B
in 25 minutes**

Eluent A: 20 mM KH_2PO_4 pH 4.5 *or* 7.0

Eluent B: 85 % (v/v) acetonitrile, 15 % (v/v)
10 mM KH_2PO_4 pH 4.5 *or* 7.0

Tryptic digest separation using ZIC™-HILIC



Retention orthogonal to a reversed phase separation

Linear gradient 0-50 % of Eluent B in 40 minutes


Eluent A: 81 %(v/v) acetonitrile, 19 %(v/v) 15 mM KH₂PO₄ pH 4.5

Eluent B: 30 %(v/v) acetonitrile, 70 %(v/v) 20 mM KH₂PO₄ pH 4.5


Flow rate: 0.8 mL/min, UV detection 214 nm
Injection volume: 20 µL

ZIC™-HILIC column (4.6×150 mm, 5 µm)


Effect of ionic strength on cationic peptides



100 mM



50 mM



30 mM

Peptides

(1) Asp-Arg-Val-Tyr-Ile-His-Pro-Phe; (2) Arg-Pro-Pro-Gly-Phe-Ser-Pro-Phe-Arg; (3) Gly-His-Lys.



Linear gradient ranging from 10-70 % eluent B in 20 minutes

Eluent A: 85 % (v/v) ACN, 15 % (v/v) 10 mM KH_2PO_4 pH 4.5.

Eluent B: 100 mM or 50 mM or 30 mM KH_2PO_4 pH 4.5;

ZIC™-HILIC column (4.6×150 mm, 200Å, 5 μm)

ZIC™-HILIC separation of plant hormones



Linear gradient 0-15 % in 10 minutes and 15-27 % in 30 minutes of Eluent B


Eluent A: 100 % (v/v) acetonitrile

Eluent B: 6.5 mM ammonium acetate pH 5.5

ZIC™-HILIC column (3.0×150 mm, 200Å, 5 μm)

With kind permission from Anders Nordström (Dept. of Forest Genetics and Plant Physiology, SLU, Umeå, Sweden)

ZIC™-HILIC separation of a plant extract



General structure of
sapogenin glycoside, i.e.,
steroidal saponin


R = sugar moiety; glucose,
galactose, pentose etc.

Linear gradient 10-100 % in 120 minutes

Eluent A: 96 % (v/v) acetonitrile 4 % (v/v)
10 mM ammonium acetate pH 5.6

Eluent B: 75 % (v/v) acetonitrile 4 % (v/v)
10 mM ammonium acetate pH 5.6

ZIC™ separation of proteins



- ❖ Loading capacity of lysozyme ≈ 42 mg/mL column volume
- ❖ Recovery of lysozyme $>95\%$


**Gradient: 30 % of Eluent B for 3 minutes
30-100 % of Eluent B in 4 minutes**

Eluent A: 5 mM phosphate buffer pH 6
Eluent B: 5 mM phosphate buffer pH 6
250 mM NaCl


ZIC™-Si column (4.6x50 mm, 200Å, 10 μ m)

Generally milder conditions than for SCX columns

Stability of the ZIC™ stationary phase



pH 2 at 50 °C



pH 8 at 50 °C

Conclusions

- ❖ The ZIC™-HILIC column is suitable for peptides and polar compounds poorly retained on RPLC phases.
- ❖ ZIC™ can be employed as a complement to cation-exchange chromatography of charged peptides.
- ❖ Volatile buffers can be employed.
- ❖ The electrostatic effect may be utilized to improve the selectivity and resolution in the ZIC™-HILIC mode
- ❖ ZIC™ is a useful tool, soon in everyone's hand... (*sic!*)

